---

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The electronical baby scale has been designed to follow the evolution of the child's weight and height during all phases of his growth.

## BABY SCALE

Use from birth to monitor the weight of the infant, quantifying the amount of milk assimilated to each feed reliably with a division of 1/2 oz (10g).

The TARE function allows the baby to be weighed after placing a cloth or towel on the tray. The baby scale is also equipped with a WEIGHT-BLOCK function which allows you to display a constant weight despite the fact that the child keeps moving.

As soon as the child can stand up alone, the weighing operations can be carried out without the tray up to the maximum weight of 44,40 pounds (20 kg).

## MEASURING ROD

A built-in height makes it possible to measure the child up to 1m50 high. Reassuring, the platform and feet of the scale are non-slip.

**Attention ! Read carefully the instructions and warnings before using the equipment and keep it in a safe place for the future reference.**

## TIPS FOR USE:

- Place the scale on a flat and hard surface. A soft and uneven floor (e.g. carpets, linoleum) could cause an alteration to the weight reading. The tray must be well positioned horizontally.
- The Baby Scale with tray is designed to weigh a lying infant. If the baby is weighed sitting down, then it should be positioned in the middle of the tray. Pay attention the baby stands still while the scale measures the weight.

**CAUTION** - To avoid any accident, the baby placed on the tray must remain under constant adult supervision.

- If the scale is used without tray, the child must be able to stand alone on the scale and not weigh more than 20 kg. Pay attention that the baby does not lean on anything while being measured, as this could alter the reading.

**CAUTION** - Never leave the child unattended.

## INSTALLING BATTERY:

Use a 9V alkaline battery only. Remove the battery cover on the underside of the scale. Insert the battery respecting the polarities and refit the cover. Do not throw away the battery together with your domestic rubbish. Take them to a used batteries bin.

## BUTTONS FUNCTION:

**ON BUTTON:** The "ON" button turns the scale on.

**OFF BUTTON:** The "OFF" button turns the scale off.

**UNIT BUTTON:** Press the "UNIT" button to select the unit of measure you require : for weight (kg or pounds) and for size (cm or ft.in).

This scale has a capacity of 20 kg, division 10 g. or 44.40 pounds, division 1/2 oz.

After turning the scale on, if the display shows "0.00" the weight is expressed in kg and g, if the display shows "0: 0" the weight is expressed in pounds and ounces.

## HOLD BUTTON (BLOCK WEIGHT FUNCTION):

Due to the baby's movements, the weight shown on the display changes all the time. When the baby is laid on the tray, press the "HOLD" button.

The text "HOLD" blinking on the display indicates that the scale is calculating automatically the average of the detected weights. When the text "HOLD" stops blinking on the display, the weight of the child is displayed.

## TARE BUTTON:

If the baby is laid on a towel or cover, use the TARE function.

Turn the scale on with the "ON" button.

The display shows "8888". Wait until the display shows "0.00" then place the towel or cover on the tray. Press the "TARE" button and wait until the display shows "0.00"; the scale is now ready for use.

## MODE BUTTON:

Press the "MODE" button to select the baby scale or the measuring rod function.

## USING THE SCALE WITH THE TRAY:

Remove the 4 rubber plugs. Fix the tray to the scale using the supplied locking nut (FIG A). Press the "ON" button. The display shows "8888". Wait until the display shows "0.00"; the scale is now ready for use. Lay the baby on the tray and wait until the display shows its weight. Press the "OFF" button to turn the scale off, otherwise the scale turn off automatically after 2 minutes without being used.

## USING THE SCALE WITHOUT THE TRAY:

This baby scale can also be used without the tray, when the baby can stand up on his own. Remove the locking nut and the tray (keep for future use). Replace the 4 rubber plugs. Turn the scale on and wait until the display shows "0.00". Before the zeros disappear, get the child to stand on the scale with the feet parallelly positioned and the weight must be equally distributed (FIG B).

## USE OF THE MEASURING ROD:

The scale is equipped with the height-measurement function.

When the child can stand safely alone, it is possible to measure his/her height to monitor his/her growth. Remove the tray from the scale. Switch on the scale by pressing the "ON" button. The display shows "8888". Wait until the display shows "0.00". Let the child get on the platform as described above and press the "MODE" button. On the display appears "0.0 cm". Pull upwards for 2/3 cm the knob positioned over the display and insert it in the height fixing accessory according to figure C. Pull the accessory upwards slowly and without sudden movements until it is positioned on the head of the child (figure D).

Now, the display shows the height in cm; the text "HOLD" indicates that the data is definitive (average of the perceived heights). It is possible to measure the height in feet (ft.in) by pressing the "UNIT" button before measuring the height of the child or just after the data has appeared on the display. After having ended the height measurement procedure, accompany slowly the knob.

#### MAINTENANCE:

- Before first use, clean the tray with warm soapy water, rinse thoroughly and dry. After use, clean for hygienic reasons.
- Clean the platform with a damp cloth. Clean the bottom of the scale with a soft and dry cloth.
- Do not use abrasive sponges and chemicals or cleaners.
- Take care that no liquids penetrate into the scale, as it is not waterproof.
- Do not immerse the device in water.
- Store the baby scale in a dry environment and make sure it is turned off.
- Remove the battery if the baby scale is not intended to be use for a long time.

#### PRECAUTIONS:

Take care of your scale, as it is a precision instrument. Therefore:

- Do not expose the baby scale to extreme temperatures, humidity, shock, dust, direct sunlight or near a heat source.
- Do not drop or knock the scale.
- Do not jump on the scale.
- Do not store anything on the baby scale or subject it to an extended load; it could damage the measuring mechanism.
- Do not open the equipment for any reason; opening or altering the equipment permanently invalidates the warranty.
- Do not leave or use the scale in very damp surroundings, humidity above 85% or contact with water and other liquids could cause corrosion and effect the precision of the scale.

#### SYMBOLS:

- "LO" indicates that the battery is flat. Replace it.
- "ERR" shows the weight is too much.
- The scale is ready for use when the display shows "0.00"; if the instructions have not been followed correctly, the display shows "ERR" or "----".

#### TECHNICAL DATA:

Range of measurement : 1 kg - 20 kg / 30 cm - 150 cm Division: 10 g

Battery : 1 alkaline battery 6LR61 series

Accuracy within : +/-1% + 1 division of applied weight

Working conditions : +10°C +40°C RH 85%

Storage conditions : +10°C 50°C RH 85%

Dimensions : 28.5 x 55 x 16 cm

Weight : 2.9 kg

#### WARRANTY:

This equipment is covered by a 2-years warranty starting from the purchase date.

It covers any defect resulting from a material or manufacturing which is not commercially acceptable. dBb Remond products are designed for domestic use only.

The warranty does not apply if the damage is caused by:

- accidental events, misuse, negligence or improper use of the product. And in case of anomalies caused by humidity, fire, or natural disasters;
- use of parts or accessories not recommended by the manufacturer;
- unauthorized repair or modification performed by third parties;
- unsuitable transport of the item without adequate packaging;
- in the device is not accompanied by the purchase invoice.

Warranty does not apply to parts subject to tear and wear along the use, as well as to batteries, when they are included.

In case of faults, please contact to the dealer where you purchased the device first; do not send anything directly to dBb Remond.

The manufacturer cannot be made liable for any damages that may be caused, directly or indirectly, to people, things and pets resulting from the inobservance of the provisions contained in the suitable user manual, especially with reference to the instructions related to installation, use and maintenance of the equipment.

Being constantly involved in improving its own products, dBb Remond reserves the right to modify its product, totally or partially, with reference to any production need, without that this involves any responsibility whatsoever by dBb Remond or its dealers.

#### DISPOSAL PROCEDURE (Dir. 2012/19/UE-WEEE)

The symbol indicates the separated sorting of electric and electronic equipment, to ensure that it is recycled in a environmentally safe manner.

At the end of life of the device, do not dispose it as mixed household waste, but dispose it in a specific collect for recycling electrical and electronic equipment.

S.E.B.I.R  
54 Rue Rouget de Lisle  
93697 PANTIN cedex - France  
www.dbb-remond.fr

NOTE

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

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54 Rue Rouget de Lisle  
93697 PANTIN CEDEX

+33 (0)1 49 15 96 40

[www.dBb-remond.fr](http://www.dBb-remond.fr)

[contact@dbb-remond.fr](mailto:contact@dbb-remond.fr)

[www.facebook.com/dBb.Remond](https://www.facebook.com/dBb.Remond)

[www.instagram.com/dbbremond](https://www.instagram.com/dbbremond)